

1.5 ECO



Next Day Delivery
Best Price Available
Quality Lenses


ECO 1.5
SINGLE VISION
SPHERIC

UNCOATED

STOCK

PRODUCT
CODE

155




n _e	1.50
V	58
	1.32
	355nm
UVA	93%
UVB	100%

STOCK

Ø60/Ø65/Ø70

Ø55



ECO 1.5
SINGLE VISION
SPHERIC


HARDCOAT

STOCK

NON
TINTABLE


PRODUCT
CODE

155



n _e	1.50
V	58
	1.32
	355nm
UVA	93%
UVB	100%

STOCK



ECO 1.5
SINGLE VISION
SPHERIC


HMC

STOCK

NON
TINTABLE


PRODUCT
CODE

155



n _e	1.50
V	58
	1.32
	355nm
UVA	93%
UVB	100%

STOCK



1.6 ECO



Best Available Price
Stock 2 Day Delivery
Quality Lenses

ECO 1.6
SINGLE VISION
SPHERIC


HMC

STOCK

NON
TINTABLE

PRODUCT
CODE

156



n _e	1.60
V	36
	1.36
	380nm
UVA	99%
UVB	100%

STOCK



1.6 Stock has a delivery date of 2 days

ECO 1.6
SINGLE VISION
ASPHERIC


HMC

STOCK

NON
TINTABLE

PRODUCT
CODE

157



n _e	1.60
V	36
	1.36
	380nm
UVA	99%
UVB	100%

STOCK



1.6 Stock has a delivery date of 2 days

ECO 1.6
SINGLE VISION
ASPHERIC


HMC

RX

NON
TINTABLE


PRODUCT
CODE

157




n _e	1.60
V	36
	1.36
	380nm
UVA	99%
UVB	100%

STOCK



1.6 Rx has a delivery date of 4-5 days

ECO




Next day delivery
Best price available
Quality Lenses

*Precal, special treatments and tints not available
*Net prices



ECO




2007

ESSILOR IRELAND SALES,
LOUGHMORE AVENUE,
RAHEEN BUSINESS PARK, RAHEEN, LIMERICK.
TEL: 061 228 553
FAX: 1800 776 776
e: IEorders@essilor.ie
WWW.ESSILOR.IE
WWW.VARILUX.IE

M15ECO

ECO



Essilor
ECO Lens Range 2007

