


Ormix™ 1.6

The material you can rely on



Material features

	PRODUCT FEATURES
INDEX	1.6
ABBÉ NUMBER	41
DENSITY	1.30
UV CUT OFF	400nm

Ormix 1.6 is a material you can recommend to wearers with confidence. Ormix 1.6 lenses perfectly meet wearers' needs: they are discreet and are suitable for any occasion.

Thin lens

Thanks to its refraction index, Ormix 1.6 is up to 30% thinner than Orma (average thickness of a single-vision lens edged to 50mm)

UV Protection

Ormix 1.6 lenses offer 100% protection against harmful UVA & UVB* rays: 400nm cut off

*UV rays can damage the retina just as they can harm the skin. Even non-tinted Ormix 1.6 lenses ensure a perfect UV cut off.

Ultimate optical correction

Ormix 1.6 material offer excellent contrast and clear images (V=41).

Suitable for the most popular frames

Ormix 1.6 gives you perfect mounting every time – including rimless frames
- Easy to edge
- Excellent resistance to drilling

LENS ORDERING BY TELEPHONE
061 228553

LENS ORDERING BY FREE FAX
1800 776 776

LENS ORDERING BY EMAIL
IEORDERS@ESSILOR.IE



A wide range

Single Vision Ormix 1.6 spherical and aspherical lenses are available in prescription and in stock, in the most extensive power range. The optical design has been enhanced to create thinner lenses and improve quality.

Ormix 1.6 lenses are also available in all the latest Varilux design: Varilux Ipséo and Varilux Physio.

Varilux Ellipse and Varilux Comfort are available too.

Both Single Vision and Varifocals lenses come as clear or Transitions.

Ormix 1.6 lenses are available with Essilor' state of the art coating Crizal Alizé with Scotchgard Protector, Crizal Alizé and Supra. See table below for full availability.

Ormix 1.6 lenses can be tinted with Essilor tints.

Crizal' Alizé™ with



	CRIZAL ALIZÉ WITH SCOTCHGARD	CRIZAL ALIZÉ	SUPRA
ORMIX SINGLE VISION			
SPHERIC CLEAR STOCK		•	•
SPHERIC CLEAR RX	•	•	•
ASPHERIC TRANSITIONS RX	•	•	
ASPHERIC CLEAR STOCK		•	
ASPHERIC CLEAR RX	•	•	
ORMIX VARIFOCALS			
VARILUX IPSEO CLEAR & TRANSITIONS	•		
VARILUX PHYSIO CLEAR & TRANSITIONS	•	•	
VARILUX ELLIPSE CLEAR & TRANSITIONS	•	•	
VARILUX COMFORT CLEAR & TRANSITIONS	•	•	